

Far Infrared Heating Pad

The Far Infrared Heating Pad is used for superior Far Infrared (FIR) Technology- a proven form of heat therapy that works more efficiently and effectively than traditional heating techniques.

Best for pain/stress relief-muscle relaxation , Improvement of blood circulation and promote healing.

Far Infrared

- Far infrared is part of the spectrum of light emanating from the sun.
- Far infrared is not visible to the eye, but you feel it as heat penetrating into your skin.
- Far infrared emits safe far-infrared rays that penetrate deep into the joint and surrounding tissues which increase blood flow, reduces inflammation and decreases pain.
- FIR was called as *Light Life*.



3-Stage Temperature Cycle

12V, 1.5A, 30 minute auto shut-off for safety

- Low Temperature (40-45°C)
- Medium Temperature (45-50°C)
- High Temperature (50-55°C)

Features:

- Comes complete with different braces, Temperature Controller, Heating Pad, AC/DC Adaptor.
- There are Shoulder, Elbow, Hand & Wrist, Back, Thigh, Calf, Knee, Ankle, Foot pads, which is ergonomically designed to fit the body part.
- Designed from 100% carbon fiber fabric which in itself is an electrically conductive material that naturally emits infrared energy to support therapeutic application.
- Ergonomic, breathable design provides support, comfort and 100% coverage from a second injury.

Instructions:

- Remove heating pad before washing, hand wash the brace by neutral cleaning solution and hang dry.
- Do not bleach, dry clean, iron, wring or tumble dry.
- For the back/hip/rib brace, the USB port can be supplied, heating therapy by rechargeable battery more convenient.



Shoulder



Wrist



Back/ Hip/ Rib



Thigh



Knee



Ankle



Foot